
SafeSpeed Warming Media kit

Open System

COD.: 38090

ESPAÑOL

USO PREVISTO

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ENGLISH

INTENDED USE

SafeSpeed Warming Media Open System is designed for ultra-rapid warming of human oocytes and embryos previously cooled by vitrification.

PRODUCTS FOR WARMING INCLUDED IN THE KIT:

4 x Vial 1 (2ml): Thawing Solution (TS)
1 x Vial 2 (2ml): Diluent Solution (DS)
1 x Vial 3 (2ml): Washing Solution (WS)

OTHER MATERIALS NEEDED FOR WARMING NOT INCLUDED IN THE KIT:

- Vitrification straws (with the samples)
- Metal clamps
- Sterile Petri dish, suitable for embryology use
- Liquid Nitrogen container (15 cm depth)
- Aspiration system (Stripper Pipette) *
- Microscope
- Stopwatch or timer
- Micropipette (100-1000 µL)
- Straw cutter: sharp scissors

***Note:** Use a pipette with a suitable internal diameter for oocytes or embryos. 170 µm capillaries are recommended for oocyte vitrification and 170 µm capillaries for embryos. For early blastocysts (up to type 3 according to Gardner or ASEBIR classification), 220 µm capillaries are recommended and 290 µm capillaries for expanded blastocysts. This allows optimizing the volume of the solutions to achieve the best dilution condition and obtain a high survival rate.

CHEMICAL COMPOSITION

Base components (Vial 1: Washing Solution – WS)

Phenol Red	Na-Pyr	NaOH	Sucrose
PBS	HPC	HCl	

Vial 1: Thawing Solution (TS)	Vial 2: Diluent Solution (DS)
Sucrose 1 M	Sucrose 0.5 M

QUALITY CONTROL TESTING

Each solution of the SafeSpeed Warming Media kit Open System is sterilized by filtration and confirm a sterility assurance level (SAL) 10⁻⁶. Each lot of SafeSpeed Warming Media Open System is tested for:

- Endotoxins by LAL methodology ≤ 0.5 EU/ml
- Mouse Embryo Assay (one-cell) ≥ 80 % expanded blastocyst after 96 h
- Sterility test according to UNE EN ISO 11737-2
- pH test: 7.2-7.6

All results are reported on a lot-specific Certificate of Analysis which is available upon request.

STORAGE AND STABILITY INSTRUCTIONS

The products are aseptically processed and supplied sterile. The solutions should be stored in their original, unopened vial and refrigerated at a temperature between 2-8 °C protected from (sun) light. Do not use media that have not been stored under these conditions. When stored as directed by the manufacturer the product is stable until the expiry date shown on the vial label. The product is supplied in single-use vials.

Do not freeze before use.

WARNINGS AND PRECAUTIONS

- Read the instructions for use before using the product.
- Do not re-sterilize.
- Do not re-use. This product is for single use and not to be reused due to the risk of contamination and infection.
- Not injectable.
- Do not use any vial if shows evidence of damage, leaking, suspended particles, turbidity, or has changed colour. Discard the product in accordance with the applicable standards.
- The product should not be used if the sterile package is found to be incomplete or open.
- Do not use and discard the product if it has passed its expiry date.
- Dispose of excess (not used) after warm-up.
- Each media is for same patient use.
- This product is intended to be used by professional trained personal in assisted reproductive techniques that includes the intended application for this product.
- To avoid contamination problems, aseptic handling techniques should be used.
- Additional precaution, it is recommended to carefully examine the vitrification straw when removing it from its package and check it does not show cracks or ruptures.
- Dispose the product as residual and non-reusable solutions.
- Dispose the containers and packaging as indicated to containers with solutions inside.
- It is recommended to use sealed vitrification instruments.
- Currently, clinical research literature indicates that the long-term effects of vitrification on oocytes and embryos, and on new-born infants born after cryopreservation remain unknown.
- The User facility of this device is responsible for maintaining product traceability and must comply with national regulations regarding traceability, where applicable.
- It is strongly recommended that the patient's name and the batch number of the SafeSpeed Media product being administered be recorded to ensure the patient-product relationship. Vecmedical Spain S.L. will keep the information regarding the product batch for 20 years.

CONTRAINDICATIONS

This product contains Ethylene glycol (EG) and dimethyl sulfoxide (DMSO), the user must take the necessary precautions to ensure that the patient is not sensitive to any of these components.

I INSTRUCTIONS FOR USE

Make sure that all media have their content homogenized before using them by gently inverting them twice. It is recommended to read all the steps of the vitrification procedure before starting the process. The process should be performed at room temperature (~22-27 °C).

PREPARATION FOR WARMING

1. Heat the TS flask (sealed) and a double-well Petri dish, or similar, (plate C) in an incubator at 37 °C (>= 1 hour) or in a CO2-free oven (do not heat in a heated plate).
2. Expose vials with DS and WS solutions to room temperature at least one hour before use.
3. Label DS, WS1 and WS2 on the petri dish cover (dish D).
4. Gently invert each vial of DS and WS twice to help homogenize their contents.
5. Using a micropipette, form the following three drops on the dish D: one drop of 100 µL from DS next to the DS mark, one drop of 100 µL from WS by the WS1 mark and one last drop of 100 µL from WS by the WS2 mark (See illustration 1). Place the dish D on the microscope and cover it.
6. Remove the Petri dish C and TS vial from the incubator. Place the Petri dish C on the microscope, gently invert the TS vial twice to help homogenize its contents and pour TS into the Petri dish C (See illustration 2) (pour enough TS to fully cover the Petri dish).
7. Adjust the focus of the microscope to the Petri dish C with the TS. Fill 90 % of the container with liquid nitrogen.

WARNING

- CAUTION:** Use metal clamps to avoid direct user contact with the liquid nitrogen.
8. Remove the vitrification device with the samples to be warmed from the storage tank and transfer it to the container with liquid nitrogen. **CAUTION:** In this process it is essential that the samples are always covered by liquid nitrogen.
 9. Set the timer. Use the timer to control the next steps that require it.
 10. Aseptically, select the vitrification straw with the samples and prepare it according to the instructions for use of the device.

CAUTION: Make sure that the capillary with the samples is constantly immersed in liquid nitrogen during handling.

11. Quickly (< 1 second), remove the vitrification straw with the samples and immerse it in the TS medium in Petri dish C.
12. Follow the directions of the specific rewarming protocol of the cryopreservation device you are using to release the samples inside the TS medium
13. Leave the samples in the TS medium for a maximum of 1 minute (see Figure 3). If possible, perform several washes to remove any residual VS with a micropipette of appropriate diameter, without touching or inserting it into the micropipette.

DILUTION

14. Aspirate with the pipette tip the oocyte/embryo with a small volume of TS and gently place it at the bottom of the drop of 100 µL from DS of the dish D. Leave it for 3 minutes (See illustration 4a-b).

WASHING

15. Aspirate the oocyte/embryo from the DS with a small volume of DS, and gently place it at the bottom of the 100 µL drop of WS1. Leave it for 5 minutes (See illustration 5a-b).
16. Aspirate the oocyte/embryo with the minimum volume of WS1 with the pipette tip and transfer it onto the drop surface of 100 µL of WS2 (See illustration 6a-b).
17. When the oocyte/embryo falls to the bottom of the WS2 drop, repeat the work on the WS2 drop. For one minute (See illustration 6c).
18. Transfer the oocyte/embryo to a culture plate with the appropriate culture medium. Incubate the oocyte/embryo in a CO2 incubator at 37 °C for complete recovery before microinjection/transfer. **Note:** a two-hour recovery of the oocyte/embryo is recommended.

SAFETY INFORMATION

The Safety Data Sheet is available to the user. Please contact the manufacturer for the MSDS.

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Vecmedical Spain SL.
Can Milans Nave 9, Montcada i Reixac,
08110 Spain.
DISTRIBUTED BY: SafePreservation S.L.

REVIEW 06/2023