

DAKEWE

HistoFlo™ Embedding Cassettes



Dakewe's HistoFlo™ cassettes are engineered to optimize print quality when used with laser cassette printers from Dakewe. The cassettes allow efficient processing and specimen safety. HistoFlo™ cassettes are made of P.O.M materials and resistant to common processing chemicals and heat.



Pore: Slotted
Open: Front open

002019□□ 002009□□
90 pcs/tube 75 pcs/tube
loose packed attached lids
lids



Pore: Slotted
Open: Back open

002119□□ 002109□□
90 pcs/tube 75 pcs/tube
loose packed attached lids
lids




Pore: Square
Open: Front open

012019□□ 012009□□
90 pcs/tube 75 pcs/tube
loose packed attached lids
lids



Pore: Slotted
Open: Front open

001009□□
One-piece
200 cassettes/tape



Pore: Square
Open: Front open

011009□□
One-piece
200 cassettes/tape



Pore: Mesh
Open: Front open

031000□□
One-piece
200 cassettes/tape



Pore: Round
Open: Back open

022129□□
90 pcs/tube
no lid



Pore: Mesh
Open: Back open

032109□□
80 pcs/tube
attached lids

Simply select your cassette and configuration number above, add the two digit color code from the right.
Example: Green Mesh Biopsy Cassettes = 032109GN



Pink
PK



White
WH



Yellow
YW



Green
GN