


NEW ROTOR Click Spin Lid (CSL)




Finally !
The New High security sealed Rotor

The only Biological containment rotor that cannot be used **without its lid**



Safer for users and hazardous samples

- No rotor lid = Incorrect operation
- If the Rotor lid is missing:
the AFI centrifuge will automatically not operate
- A warning message will appear on the screen
- Rotor lid = Correct operation
=> **Biological safety guarantee**



Second element to ensure a secure lid closure

- You have a visual indicator to secure the biological containment
- Red = lid is not secured
- Turn until green alignment is reached
- The lid is now secured

